

Acceleration levels h/p/cosmos running machines



acceleration level (nr.)	0 to maximum speed (s)	acceleration in m/s ²	maximum treadmill speed in km/h								
			10	22	25	30	40	45	50	60	80
1	131		0,021	0,047	0,053	0,064	0,085	0,095	0,106	0,127	0,170
2	66		0,042	0,093	0,105	0,126	0,168	0,189	0,210	0,253	0,337
3	33		0,084	0,185	0,210	0,253	0,337	0,379	0,421	0,505	0,673
4	16		0,174	0,382	0,434	0,521	0,694	0,781	0,868	1,042	1,389
5	8		0,347	0,764	0,868	1,042	1,389	1,563	1,736	2,083	2,778
6	5		0,556	1,222	1,389	1,667	2,222	2,500	2,778	3,333	4,444
7	3		0,926	2,037	2,315	2,778	3,704	4,167	4,630	5,556	7,407

