



There is an option to combine the h/p/cosmos arm support to a detachable additional keyboard (6 buttons). Start, stop, speed and elevation are adjustable for the subject from a safe and ergonomic position.

The additional keyboard can be taken from its mounting. With the long connection cable the therapist has the possibility to control the running machine, e.g. from a kneeling position.

 [Additional stop button](#)

 [h/p/cosmos arm support](#)

 [Retro-fitting](#)

 [Printout](#) [*.pdf]

Mounts for the h/p/cosmos additional keyboard



Two mounts for the h/p/cosmos additional keyboard are available. One for the handrail and a second one, to fix the remote control beside the motor hood.

order numbers	
Additional keyboard incl. 2 m cable	cos10106
Mount for additional keyboard arm support	cos10111
Mount for additional keyboard handrail	cos14135
Mount for additional keyboard motor left	cos14327
Mount for additional keyboard motor right	cos11750
Additional cable 2 m	cos12922