




The handrails adjustable in height and width offer best possible safety at best attainable freedom of body motion. The adjustable handrails are useful especially for exercise and all therapies with children and smaller people and to insecure or elderly people.

-  [Joints with scales](#)
-  [Stabilising elements](#)
-  [Adjustment](#)

Optionally available with extra charge for the running machine series ...

-  [h/p/cosmos mercury](#)
-  [Printout \[*.pdf \]](#)

The handrails are to be adjusted individually at the two [joints with scales](#) at the back end. The two [stabilising elements](#) at the front are responsible for the fixation of the adjusted position.




If using the adjustable handrails, the user terminal is placed at the crossbar in front of the testing person.



Because of the tremendous variability and broad adjustments the handrails are applicable for most body heights.

Especially for children they offer additional safety in therapy and exercise - but: the [h/p/cosmos safety arch](#) with fallstop must be used as well.



| | |
|--------------------------------|----------|
| order number | |
| h/p/cosmos handrail adjustable | cos10030 |