



The model series [h/p/cosmos mercury 4.0](#) can optionally be equipped with short handrails. They offer access for the therapist or physician to the patient and best freedom of body motion during training.

Optionally available with extra charge for the running machine series ...

 [h/p/cosmos mercury](#)

 [Printout](#) [ \*.pdf ]




The short handrails make the running machine extremely compact. Therefore it is ideal even for smaller surgeries, e.g. [h/p/cosmos mercury lt med 4.0](#) for cardiology. In this application the control of the running machine usually occurs via interface from the ECG system. In this case a user terminal is not necessary.



Athletes are enthusiastic about the extrem freedom of body motion when training with [h/p/cosmos mercury 4.0](#) with short handrails. For motion analysis the short handrails always offering a wide field of view.

At the [h/p/cosmos mercury 4.0](#) with short handrails, the user terminal is not placed on the right hand side handrail. For technical reasons it is placed at the crossbar in the front.



order number	
<a href="#">h/p/cosmos handrail short 150 [mercury]</a>	cos11456